

**AddOn®**

For the correction of pseudophakic eyes


# AddOn<sup>®</sup> sulcus lens

The AddOn is the IOL platform for vision enhancement and correction of pseudophakic eyes. The patented sulcus lens constitutes a milestone in functional patient care.

The AddOn lens is implanted in the sulcus as an addition to the Basis IOL fixed in the capsular bag, and is compatible with commonly used capsular bag IOLs, regardless of design or material.

## Basis IOL

in the capsular bag




## AddOn<sup>®</sup> IOL

in the sulcus

# AddOn<sup>®</sup> Areas of Application


## Optimizing

refraction results in pseudophakic eyes. The full optical correction spectrum is available and can be individually adapted to the patient.




## Upgrading

to multifocal (trifocal) or EDOF function in case of a subsequent patient request for more spectacle independence.



## Planning spectacle independence


with the Liberty<sup>2</sup> dual-lens procedure. Spectacle independence becomes plannable and predictable with a Basis IOL in the capsular bag and a trifocal AddOn in the sulcus.




## Optical range

	Distance	Near	Intermediate	Imaging
1  AddOn <sup>®</sup> spheric	●			
2  AddOn <sup>®</sup> toric	●			●
3  AddOn <sup>®</sup> EDOF	●		●	
4  AddOn <sup>®</sup> trifocal	●	●	●	
5  AddOn <sup>®</sup> toric-trifocal	●	●	●	●


## The AddOn<sup>®</sup> features




**Sulcus fit**  
through 4 flexible haptics



**Rotational stability**  
through non-torque design



**No IOL touch & cell stop**  
through convex-concave optic



**No iris capture**  
through square design

## The AddOn<sup>®</sup> Benefits

The AddOn offers significant advantages for individual and safe patient care.

- **Patented sulcus design<sup>(1)</sup> for safe and stable positioning in the ciliary sulcus<sup>(3)</sup>**
- **Full spectrum of refractive corrections<sup>(2)(3)(9-12)</sup>, in particular trifocal upgrades<sup>(4)</sup>**
- **Safe procedure<sup>(2)(5)(6)</sup> based on a proven technology**
- **No compromise of optical quality through the additional AddOn interfaces<sup>(9)</sup>**
- **Non-traumatic reversibility<sup>(9)(10)</sup>**
- **Implantable irrespective of lens status of the pseudophakic eye<sup>(5)(7)</sup>**
- **Option of implantation in one procedure or two subsequent procedures<sup>(9-12)</sup>**
- **Suitable for eyes with silicone oil tamponade or previous keratoplasty<sup>(8)</sup>**

### References

- (1) Patent #EP 12 181 744 (2013)
- (2) Gundersen KG, Potvin R. A review of results after implantation of a secondary intraocular lens to correct residual refractive error after cataract surgery. *Clinical Ophthalmology* 2017;11, pages 1791-1796.
- (3) Gundersen KG, Potvin R. Refractive and visual outcomes after implantation of a secondary toric sulcus intraocular lens. *Clinical Ophthalmology* 2020; 14:1337-1342.
- (4) Palomino-Bautista C (et al). Spectacle independence for pseudophakic patients - Experience with a trifocal supplementary add-on intraocular lens. *Clinical Ophthalmology* 2020; 14:1043,1054. eCollection 2017.
- (5) Mrochen M (et al). Sicherheit moderner Zweilinsensysteme. Modellbasierte Untersuchung hinsichtlich Lagetoleranz und Auswirkungen auf das optische System. *Ophthalmologie (Sonderveröffentlichungen, Teile 1-3)* 2013.
- (6) Srinivasan S. Implantation of Scharioth macula lens in patients with age-related macular degeneration: results of a prospective European multicentre clinical trial. *BMJ Open Ophthalmol.* 2019 Jul; 4(1):e000322.
- (7) Reiter N, Werner L, Guan J, Li J, Tsaousis KT, Mamalis N, Srinivasan S. Assessment of a new hydrophilic acrylic supplementary IOL for sulcus fixation in pseudophakic cadaver eyes. *Eye (Lond)*. 2017 May;31(5):802-809.
- (8) Hassenstein A, Niemeck F, Giannakakis K, Klemm M. Torische AddOn-Intraokularlinsen zur Korrektur hoher Astigmatismen nach pseudophaker Keratoplastik. *Der Ophthalmologe.* 2017 Jun;114(6):549-555
- (9) Khoramnia R, Baur ID, Yan W, Łabuz G, Auffarth GU Comparison of a Presbyopia-Correcting Supplementary Intraocular Lens Combination and a Capsular-Bag Lens: An In Vitro Study. *Diagnostics* 2023, 13, 1482. <https://doi.org/10.3390/diagnostics13081482>
- (10) Gundersen KG, Gjerdrum B, Potvin R. Efficacy of a Secondary Trifocal Sulcus IOL in Providing Near and Intermediate Vision in Patients with Prior Myopic Laser Vision Correction and Cataract Surgery. *Clinical Ophthalmology* 2022;16:2219-2226.
- (11) Albayrak S. et al. Visual performance and patient satisfaction following the implantation of a novel trifocal supplementary intraocular lens. *European Journal of Ophthalmology* 2021;31(5):2346-2352.
- (12) Harrisberg BP, Chua AW, Chua MJ, Taher A Comparison of Primary Duet Lens Procedures: In-The-Bag Monofocal with Sulcus Multifocal, and Standard Single Multifocal Lens for Cataract Surgery. *Clinical Ophthalmology* 2023;17:273-282. <https://doi.org/10.2147/OPTH.S396472>

## Technical Data

	 AddOn® spheric	 AddOn® toric	 AddOn® EDOF	 AddOn® trifocal	 AddOn® toric-trifocal
<b>Sphere/SEQ</b>	-10.0 D to +10.0 D	-10.0 D to +10.0 D	-5.0 D to +5.0 D	-5.0 D to +5.0 D	-3.0 D to +3.0 D
<b>Cylinder</b>		+1.0 D to +11.0 D			+1 D to +4.5 D
<b>Addition</b>			+1.75 D	+2.25 D; +3.00 D	+3.0 D

<b>Type and material</b>	One-piece sulcus lens made of biocompatible hydrophilic acrylic (water content 25 %), UV absorber
<b>Refractive index</b>	1.46
<b>Abbe number</b>	58
<b>Dimensions</b>	Total diameter: 13 mm   Optic diameter: 6 mm
<b>Optic</b>	spheric optic   convex-concave optic design
<b>Haptic</b>	4 flex-haptics   0° angulation   4 closed loops   sulcus fit
<b>Sterilization</b>	Steam
<b>Recommended injection system</b>	1stInject 2.0HB   1stInject 2.4HB
<b>Manufacturer</b>	1stQ. AddOn® is a registered trademark of 1stQ.

For use by healthcare professionals only. Please carefully read the safety information in the Electronical Instructions for Use before use (eIFU / [www.1stq.eu/eifu](http://www.1stq.eu/eifu)). In case of questions, please contact the 1stQ contact person responsible for you.

## What users around the world are saying



*"Based on the reliable stability of the 1stQ AddOn Trifocal lens, I use it often in clinical practice for increasing my options for refractive upgrades or correction."*

**Sathish Srinivasan, MD, FRCSEd, FRCOphth, FACS, United Kingdom**

*"I feel very comfortable and safe with the planned two-stage trifocality with a Basis IOL in the capsular bag and a trifocal AddOn in the sulcus. The superindication for this technique is RLE or post-LASIK cataract surgery. Refractive surprises should thus be entirely eliminated."*

**Kjell Gunnar Gundersen, MD, PhD, Norway**

*"The 1stQ AddOn IOL allows to bring patients into their visual comfort zone."*

**Erik L. Mertens, MD, FEBO, Belgium**

*"The AddOn with trifocal array is just great. In combination with the capsular bag lens it provides a high degree of adaptability to refractive situations. Handling and safety are convincing."*

**Prof. Dr. Gangolf Sauder, MD, Germany**

*"The main advantage of implanting a secondary lens with trifocal performance is its ease of reversibility at any time, [...]. This is a smart option for the patient who is particularly unsure about dealing with aberrations from multifocal IOLs."*

**Brian Harrisberg, MD, MBBCh, FRACS, FRANZCO, Australia**

# Leading begins together

**Want to find out more about the AddOn®? Or discuss  
1stQ's versatile IOL solutions?**

We are ready and happy to assist with expertise or new ideas.  
Let us drive quality, standards and the best possible individual  
vision together.

**Further questions?**

[export@1stq.eu](mailto:export@1stq.eu)

Service Box  
AddOn®



▶ Handling video



↓ AddOn website