


AddOn® SML

hydrophilic



AMD

Vision functions	IOL functions
Ultra-near function	Magnifier function
Protective function	Blue filter function (optional)
Distant function	Basis function (on request)

Features

One piece, foldable intraocular lens (IOL), biocompatible hydrophilic acrylic (25% water content), as additional IOL in combination with a Basis IOL.

- Convex/Concave Optic
- 4 Flex-Haptics
- Non Torque Design
- Square Design
- Addition
- ▶ Cell Stop/Non IOL Touch
- ▶ Sulcus Fit
- ▶ Rotational Stability
- ▶ Non Iris Capture
- ▶ +10 D; Ø 1.5 mm

AddOn® SML	AddOn® SML yellow	Sph.	Addition
A4MW00	A4MWY0	on request	+10 D