


# Basis V

## hydrophobic



### Q functional

ONV functions\*

IOL functions

Distant function

Basis function

Protective function

Blue filter function (optional)

Mesopic function

Aspherical function

### Features

One piece, foldable intraocular lens (IOL), hydrophobic acrylic. Also available as **preloaded**.

360° Cell Barrier

▶ Optimized Protection against PCO

UV Filter

▶ Retina Protection

A-constant: 118.7

Haigis-constants:  $a_0 = 1.32$ ;  $a_1 = 0.400$ ;  $a_2 = 0.100$

Manufacturer: GEMMA Medical AG, Switzerland, CE 1250

### Optical power:



Basis V clear	Basis V yellow	Basis V yellow preloaded	Sph.	Steps
B5AD00	B5ADY0		0.0 D to +9.0 D	1.0 D
B5AD00	B5ADY0	B5ABY0	+10.0 D to +30.0 D	0.5 D

\*ONV functions: The five functions of the human eye which are required to achieve optimal natural vision (ONV) (see rear page)