


AddOn®

hydrophilic • spherical



Q functional

ONV function*

IOL function

Distant function

Basis function

Features

One piece, foldable intraocular lens (IOL), biocompatible hydrophilic acrylic (26% water content), as additional IOL in combination with a Basis IOL.

- Convex/Concave Optic
- 4 Flex-Haptics
- Non Torque Design
- Square Design
- ▶ Cell Stop/Non IOL Touch
- ▶ Sulcus Fit
- ▶ Rotational Stability
- ▶ Non Iris Capture

Optical power:



AddOn®	Sph.	Steps
A4SW00	-10.0 D to +10.0 D	0.25 D

*ONV functions: The five functions of the human eye which are required to achieve optimal natural vision (ONV) (see rear page)