


AddOn[®] hydrophilic · spheric

Technical specifications

Refractive Index	1.46	Lens type	One piece sulcus lens
Abbe No.	58	Material	Biocompatible hydrophilic acrylic (25%)
		Haptic	4 Flex-Haptic
		Optic	spheric
		Optic form	convex-concave
		Filter	UV filter
		Sterilisation	Steam
		Injection system	1stInject 2.0HB, 1stInject 2.4HB



■ Basis function

Optical Power



AddOn[®] hydrophilic · spheric

Article No.	Powers	Steps
A4SW00	-10.0 D to +10.0 D	0.25 D



For the latest version,
scan the QR code.
1stq.de/addon1054

**Do you have any
questions?**
export@1stq.eu



1stQ GmbH
Konrad-Zuse-Ring 23
68163 Mannheim, Germany
www.1stq.de/en