


AddOn[®] hydrophilic · trifocal


Technical specifications

| | | | |
|------------------|------|------------------|---|
| Refractive Index | 1.46 | Lens type | One piece sulcus lens |
| Abbe No. | 58 | Material | Biocompatible hydrophilic acrylic (25%) |
| | | Haptic | 4 Flex-Haptic |
| | | Optic | spheric |
| | | Optic form | convex-concave |
| | | Filter | UV filter |
| | | Sterilisation | Steam |
| | | Injection system | 1stInject 2.0HB, 1stInject 2.4HB |




Basis function

Optical Power




Trifocal function

Addition:



AddOn[®] hydrophilic · trifocal

| Article No. | Powers | Steps | Addition |
|-------------|------------------|--------|----------|
| A4DW0K | 0.0 D | | +2.25 D |
| A4DW0M | 0.0 D | | +3.00 D |
| A4EW0M | -5.0 D to -0.5 D | 0.25 D | +3.00 D |
| A4EW0M | +0.5 D to +5.0 D | 0.25 D | +3.00 D |



For the latest version,
scan the QR code.
1stq.de/addon1051

**Do you have any
questions?**
export@1stq.eu



1stQ GmbH
Konrad-Zuse-Ring 23
68163 Mannheim, Germany
www.1stq.de/en