



Basis V

Technische Spezifikationen



hydrophob non-preloaded preloaded

A-Konstante*	Holladay 1*
118,895 (SRK/T)	1,724 (SF)
Haigis Konstanten*	Holladay 2
a0 = -0,7633	5,49 (ACD)
a1 = 0,1391	Barrett*
a2 = 0,2244	1,83 (LF); N/A (DF)
Hoffer Q*	
5,533 (pACD)	

*Stand IOL Con 03.2023

Refraktiver Index	Abbe-Zahl
1,5	50

Linsentyp	Sterilisation
1-teilige Hinterkammerlinse	EO
Material	Empf. Injektionssystem
biokompatibles hydrophobes Acryl (<0,5%)	1stInject 2.0HB, 1stInject 2.4HB / Medicel Accuject
Haptik	Filter
C-Haptik, 0° Anwinkelung	UV-Filter, Blaulichtfilter (optional)
Optik	
asphärische Optik, 360° Zellbarriere	

Lieferbereich

Basis V	Artikelnummer klar	Artikelnummer gelb	Sphäre	Schritte
non-preloaded	B5AD00	B5ADY0	0,0 D bis +9,0 D	1,0 D
non-preloaded	B5AD00	B5ADY0	+10,0 D bis +30,0 D	0,5 D
preloaded	B5AB00	B5ABY0	0,0 D bis +9,0 D	1,0 D
preloaded	B5AB00	B5ABY0	+10,0 D bis +30,0 D	0,5 D

Immer die aktuelle Version. QR-Code scannen.
www.1stq.de/basisv1001

Haben Sie noch Fragen?
+49 621 7176330 • info@1stq.de

