


611HPS Basis K

Technische Spezifikationen

hydrophil



A-Konstante*	Holladay 1*
118,3 (SRK/T)	1,39 (SF)
Haigis Konstanten*	Holladay 2
a0 = 0,95	5,14 (ACD)
a1 = 0,40	Barrett*
a2 = 0,10	1,52 (LF); N/A (DF)
Hoffer Q*	
5,15 (pACD)	

*Stand IOL Con 03.2023

Refraktiver Index	Abbe-Zahl
1,46	58

Linsentyp	Sterilisation
1-teilige Hinterkammerlinse	Dampf
Material	Injektionssystem
biokompatibles hydrophiles Acryl (25%)	1stInject 2.0HB, 1stInject 2.4HB
Haptik	Filter
C-Haptik, 9° Anwinkelung	UV-Filter
Optik	
sphärische Optik, 360° Zellbarriere	

Lieferbereich

Artikelnummer	Sphäre	Schritte
611HPS	0,0 D bis +9,0 D	1,0 D
611HPS	+10,0 D bis +30,0 D	0,5 D



Immer die aktuelle Version. QR-Code scannen.
www.1stq.de/basisk1023



Haben Sie noch Fragen?
+49 621 7176330 • info@1stq.de

