


Basis Z

Technische Spezifikationen



hydrophil

A-Konstante*	Holladay 1*
118,943 (SRK/T)	1,76 (SF)
Haigis Konstanten*	Holladay 2
a0 = 0,1657	5,52 (ACD)
a1 = 0,3113	Barrett*
a2 = 0,1631	1,85 (LF); -1,0 (DF)
Hoffer Q*	*Stand IOL Con 03.2023
5,529 (pACD)	
Refraktiver Index	Abbe-Zahl
1,46	58
Linsentyp	Sterilisation
1-teilige Hinterkammerlinse	Dampf
Material	Empf. Injektionssystem
biokompatibles hydrophiles Acryl (25%)	1stInject 1.8 bis 30,0 D, 1stInject 2.0HB, 1stInject 2.4HB
Haptik	Filter
Z-Haptik, 0° Anwinkelung	UV-Filter, Blaulichtfilter (optional)
Optik	
asphärische Optik, 360° Zellbarriere	

Lieferbereich

Artikelnummer klar	Artikelnummer gelb	Sphäre	Schritte
S1AW00	S1AWY0	-10,0 D bis -1,0 D	1,0 D
B1AW00	B1AWY0	0,0 D bis +9,0 D	1,0 D
B1AW00	B1AWY0	+10,0 D bis +30,0 D	0,5 D
S1AW00	S1AWY0	+31,0 D bis +45,0 D	1,0 D



Immer die aktuelle Version. QR-Code scannen.
www.1stq.de/basisz1014



Haben Sie noch Fragen?
[+49 621 7176330](tel:+496217176330) • info@1stq.de

