

Single piece, sterile, foldable acrylic aspheric hydrophilic sulcus intraocular lens (IOL) with UV blocker. The optic is monofocal. The IOL is non-preloaded and has to be manually loaded into a compatible injector.



CLEAR IOL

MATERIAL


Hydrophilic acrylic with 25% water content and UV blocker	
Refractive Index	1.46
Abbe Number	58

DESIGN

Optic Diameter	6 mm
Overall Diameter	13 mm
Optic Design	Aspheric - aberration neutral Convex-concave
Diffractive Zone	-
Haptic Design	Special 4 closed loops with 0° angulation
Shelf Life	5 years after sterilization

INJECTION

Estimated Incision Size	2.2 - 2.4 mm
Recommended Injector	Compatibility chart available: 1stq.eu/compatibility/



POWER RANGE

Article No.	SEQ/Powers	Steps	Cylinders	Steps	Addition
A4MW00	0.0 D*	-	-	-	+10.0 D

*Special powers on request.

Please find the online AddOn Calculator here: <https://calculators.1stq.de/AddOnCalculator>



1stQ GmbH
Konrad-Zuse-Ring 23 | 68163 Mannheim
Germany