


Single piece, sterile, foldable acrylic aspheric hydrophilic sulcus intraocular lens (IOL) with UV blocker. The optic is monofocal. The IOL is non-preloaded and has to be manually loaded into a compatible injector.



CLEAR IOL

MATERIAL

Hydrophilic acrylic with 25% water content and UV blocker

Refractive Index

1.46

Abbe Number

58

DESIGN

Optic Diameter

6 mm

Overall Diameter

13 mm

Optic Design

Aspheric - aberration neutral
Convex-concave

Diffractive Zone

-

Haptic Design

Special 4 closed loops with 0° angulation

Shelf Life

5 years after sterilization

INJECTION


Estimated Incision Size

2.2 - 2.4 mm

Recommended Injector

Compatibility chart available:
1stq.eu/compatibility/

Toric



POWER RANGE

Article No.	SEQ/Powers	Steps	Cylinders	Steps	Addition
A4TW0T	0.0 D	-	1.0 D	0.5 D	-
A4TW0T	0.0 D	-	1.5 D to 4.5 D	0.75 D	-
A4TW00	0.0 D	-	5.25 D to 8.25 D	0.75 D	-
A4TW00	0.0 D	-	9.0 D to 11.0 D	1.0 D	-
A4FW0T	-10.0 D to +10.0 D (not 0.0 D)	0.25 D	1.0 D	0.5 D	-
A4FW0T	-10.0 D to +10.0 D (not 0.0 D)	0.25 D	1.5 D to 4.5 D	0.75 D	-
A4FW00	-10.0 D to +10.0 D (not 0.0 D)	0.25 D	5.25 D to 8.25 D	0.75 D	-
A4FW00	-10.0 D to +10.0 D (not 0.0 D)	0.25 D	9.0 D to 10.0 D	1.0 D	-
A4FW00	-3.0 D to +8.0 D (not 0.0 D)	0.25 D	11.0 D	-	-

MONOFOCAL

TORIC

Please find the online AddOn Calculator here: <https://calculators.1stq.de/AddOnCalculator>



1stQ GmbH
Konrad-Zuse-Ring 23 | 68163 Mannheim
Germany